


1. Pipette **20 μ l** protease.
2. Add **200 μ l** sample.
3. Add **200 μ l** Buffer **AL**.



4. Mix by pulse-vortex for 15s and incubate at **56°C** for **10 minutes**.
5. Briefly centrifuge.
6. Add **200 μ l** 96-100% ethanol and mix by pulse-vortex for 15s.
7. Briefly centrifuge.




8. Apply the mixture to the Mini spin column.
9. Centrifuge at **8000 rpm** for **1 minute**.



10. Discard the collection tube.
11. Place the column in a new 2ml collection tube.
12. Add **500 μ l** buffer **AW1**.
13. Centrifuge at **8000 rpm** for **1 minute**.



14. Place the column in a new 2ml collection tube.
15. Add **500 μ l** buffer **AW2**.
16. Centrifuge at **14000 rpm** for **3 minutes**.
17. Place the column in a new 2ml collection tube.
18. Centrifuge at **14000 rpm** for **1 minute**.



19. Place the column in a new 1.5ml tube.
20. Add **200 μ l** buffer **AE**.
21. Incubate at room temperature for **1 minute**.
22. Centrifuge at **8000 rpm** for **1 minute**.



Add 400 μ l buffer AE to dilute the eluted DNA (3X dilution)



23. Quantify the DNA concentration.