

PHYSIOLOGY
PRACTICAL

Total and Differential Leucocytic Count **TLC & DLC**

SLIDES + HANDOUT


Red: Important

Green: Notes

Gray: Extra Information

Blue: Only boys slides

Purple: Only girls slides




437
PHYSIOLOGY TEAM



MED437
KING SAUD UNIVERSITY

Objectives :

- To be able to identify the different types of leukocytes under the microscope.
 - To practice the procedure for differential leucocyte counting.
 - To know the normal values expected for the differential white cell count.
 - To understand the use of the differential white cell count in the diagnosis of disease processes.
- 

Introduction

*Only found in boy's slides

- In the differential leucocyte count (DLC) the percentage of each type of white blood cells in the total leucocyte population is determined
- Each type of white cells, performs a different function in the battle against infections
- Each type of infection yields a different white cell picture in the blood .
- The morphology and staining characteristics of each type is peculiar and is responsible for the in specific typing.

Cell type	Erythrocyte	Lymphocyte	Neutrophil	Eosinophil	Basophil	Monocyte	Platelets
Size	6.7-7.7 µm	6-16 µm	12-14 µm	12-17 µm	14-16 µm	16-20 µm	1.5-3.5 µm
Number per litre	$3.5-5.5 \times 10^{12}$	$0.5-1 \times 10^9$	$2-2.5 \times 10^9$	$1.3-3.5 \times 10^8$	$0.5-0.8 \times 10^9$	$0.2-0.8 \times 10^9$	$150-400 \times 10^9$
Differential leucocyte count	—	20-50 %	40-75 %	1-6 %	< 1 %	2-10 %	—
Duration of development	5-7 days	1-2 days	6-9 days	6-9 days	3-7 days	2-3 days	4-5 days
Life span of mature cell	120 days	?	6 hours to a few days	8-12 days	?	Months to years	8-12 days


*only found in girl's slides

Reagents and apparatus:

- A microscope with an oil immersion objectives.
- Mineral or cedar oil
- Various dyes for staining blood films (e.g., Wright's stain and Leishman's stain) can use either stain.
- Microscope slides.

Procedure


- Prepare blood film and stain it with Wright's stain
- Examine it under the oil immersion objective lens of the microscope, and identify the different leucocytes(count about 100 cells(



WBCs Count
Between
4000-11000

This slide is VERY important


only found in girls's slides




clinical applications in the next slides are about differential count
provides clues about certain illnesses. An increase in a specific type
of Leukocytes is linked to a certain illness .

Neutrophils

40-70%
50-70%
of WBC



Erythrocytes and a neutrophil.
Wright's (Oil)




Clinical application

- **Neutrophilia:** increasing in neutrophils > acute bacterial and fungal infections (pyogenic illness)

Function:

- Active phagocytes
- Number increases
- Short term or acute infection

*It is the most common type of blood cell circulating blood




*Prominent feature: Multi-lobed nucleus

Diameter/ μm	Nucleus	Cytoplasm	Cytoplasmic Granules
<ul style="list-style-type: none"> • 14-10 	<ul style="list-style-type: none"> • Blue-violet • 6-2lobes (multilobed complex\nucleus) connected by thin strand of nucleoplasm . 	<ul style="list-style-type: none"> • Slate blue in color 	<ul style="list-style-type: none"> • Fine granules which are difficult to see closely-packed violet-pink (neutral), with various stains such as <u>Wright's stain</u>. • Small granules

Eosinophils

1 - 4/6%
1-3%
of WBC




Clinical application:

- Eosinophilia :increasing in eosinophils > Allergy and parasitic infections

Function:

- Kill parasitic worms
- Increase during allergy attacks

*It is less common in bloodstream than neutrophils



*Prominent feature: Coarse pink/red granules

Diameter/ μm	Nucleus	Cytoplasm	Cytoplasmic Granules
<ul style="list-style-type: none">• 15-10	<ul style="list-style-type: none">• 3-2 lobes• Dumbbell-shaped• often bilobed purplish• stained blue-red	<ul style="list-style-type: none">• Eosinophilic• Light pink red• Granular	<ul style="list-style-type: none">• Large• prominent ,• red• Coarse




Basophils

0-1 %
0.4- 1%

Numerous large, dark blue-violet granules that tend to be closely packed



Nucleolus is smaller round, non-segmented and stains lighter than the remainder of the cell



Clinical application:

- **Basophilia** :increasing in basophils > allergy and malignancy

Function:




- Contain histamine (vasodilator chemical) which is discharged at site of inflammation
- These large granules contain/secrete heparin and histamine
- The rarest of all white blood cells found in the blood .

*Prominent feature: Nucleus hidden by blue granules

Diameter/ μ m	Nucleus	Cytoplasm	Cytoplasmic Granules
<ul style="list-style-type: none"> • 15-10 	<ul style="list-style-type: none"> • Irregular • May be S or U shaped • Not clearly seen • Nucleus hidden behind granule (because overlaid with granules) 	<ul style="list-style-type: none"> • Basophilic • Bluish • Granular • stains dark blue 	<ul style="list-style-type: none"> • Large very coarse • Completely fill the cell • Few large blue purple granules”basophilic granules ”

Monocytes

5-10%
4 - 8%
4-6%
of WBCs




Clinical application:

Monocytosis: Monocyte count will increase with chronic infections

Function:

- Active phagocytes that become macrophages in the tissues
- Increase in number during chronic infections such as tuberculosis






*Prominent feature: Large cell with kidney shaped nucleus (no granules)

Diameter/ μm	Nucleus	Cytoplasm	Cytoplasmic Granules
<ul style="list-style-type: none"> • 20-15 • <u>The Largest WBC</u> 	<ul style="list-style-type: none"> • Single • Large • Kidney or Horseshoe shaped • dark blue purple nucleus 	<ul style="list-style-type: none"> • Abundant • Basophilic • Glass like texture 	<ul style="list-style-type: none"> • No visible granules

Lymphocytes

20-40\45 %
25 -35%
of WBC



Clinical application:

- **Lymphocytosis**(Lymphocyte count will increase with viral infections (infectious mononucleosis :

Function:

- Part of immune system
- B-Lymphocytes produce antibodies
- T-Lymphocytes involved in graft rejection, fighting tumors and viruses, and activating B-Lymphocytes.


*Prominent feature: Large nucleus taking up most of cell volume and very small cells

Diameter/ μm	Nucleus	Cytoplasm	Cytoplasmic Granules
<ul style="list-style-type: none">• Small 9-7• Some large 15-10	<ul style="list-style-type: none">• Single• Large (occupies most of cell volume (• Round• dark purple-blue nucleus	<ul style="list-style-type: none">• Thin rim (crescent) of cytoplasm• light blue color	<ul style="list-style-type: none">• No visible granules



Difference between Small and Large Lymphocytes

*only in male slides


Large Lymphocytes $\mu\text{m}15-10$



Small Lymphocytes $\mu\text{m}9-7$



DLC Equipment & Procedures : DIFFERENTIAL LEUCOCYTE COUNT



EQUIPMENTS:

1- electron microscope with an oil immersion objective 2- mineral or cedarwood oil
3- wright's stain 4- microscope slides

PROCEDURES:

A - Prepare a stained blood film with the help of Wright's stain.

“Using various dyes + microscope slide ”

B - Set the stained blood film under the oil immersion objective in an electron microscope. “Using microscope with an oil immersion objectives + Mineral or Cedar

C - Identify various types of white blood cells according to their histological characteristics.



*found in the
handout

*Values of “% of leukocytes” and “Size” may slightly vary

Summary

Blood element	%of leukocytes	Size/ μm	Cytoplasmic staining	Nucleus morphology
Erythrocyte	-	8-7	pink, no granules	none
Neutrophil	70-50	12-10	salmon-colored small granules	Segmented,- 5-2lobed
Lymphocyte	35-25	8-7	Light blue, scant amount, no granules	Single large Oval purple
Monocyte	6-4	18-16	Basophilic, no granules	Large, kidney shaped
Eosinophil	3-1	14-13	Bright red coarse granules	bilobed purplish
Basophil	1-4-0	15-14	Large, basophilic granules	Bilobed bluish black

Questions and problems:

1- What are the normal values of each different type of white blood cells?

NEUTROPHILS > 50 – 70 %

EOSINOPHILS > 1 – 3 %

BASOPHILS > 0.4 – 1 %

MONOCYTES > 4 – 6 %

LYMPHOCYTE > 25 – 35 %

2- under what conditions are the percentages of the various types of white blood cells increased?

neutrophils: will increase in acute bacterial or fungal infections.

eosinophils : will increase in parasitic infections and allergies.

basophils : will increase in allergies and malignancies.

monocytes : will increase in chronic infections.

lymphocyte: will increase in acute viral infections and malignancies.

3- What stains are used in the preparation of blood films?

[1 Leishman's stain


[2 Wright's stain

*found in the
handout


Identify the following pictures :



1




2



3



4



5

- 1Basophils
- 2Monocytes
- 3Lymphocytes
- 4Neutrophils
- 5Eosinophils

GOOD LUCK!

Team Members:
Rahaf Althnayan
Mashael Alkahtani
Dana Alrasheed

Team Leaders:
Reem ALQarni
Tareq ALOmaim

together everyone
TEAM
achieves more

References:
1- Girls' & Boys' slides
2- Handout

create
your
own
sunshine.

