

Physiology Practical For 1st Year Medical Students

(2018-2019)

Contents

List of Abbr	eviations	4
List of Figu	es	6
List of Table	es	7
Introduction	1	8
Overview o	f 1 st Year Medical Student Physiology Practical Sessions	9
Faculty and	Staff Members Involved in Physiology Practical Teaching	10
Advice for S	Students	11
Chapter 1:	Foundation Block Physiology Practical	12
Practical Sedimen	Complete Blood Count (CBC), Red blood cell indices, Erythrocyte tation Rate (ESR) and Hematocrit (Hct)	13
1.1.	Objectives	13
1.2.	Equipment	13
1.3.	Procedure	13
1.4.	Essential terminology	18
1.5.	Practice questions	19
Practical	2. White Blood Cell Count (WBC) and Differential	20
2.1.	Objectives	20
2.2.	Equipment	20
2.3.	Procedure	20
2.4.	Practice questions	22
Practical	3. Blood groups, Bleeding & Clotting Times	23
3.1.	Objectives	23
3.2.	Determination of blood groups	23
3.3.	Determination of clotting time	28
3.4.	Determination of bleeding time	30
Chapter 2:	Respiratory Block Physiology Practical	31
Practical	Simple Spirometry (Lung Volumes and Capacities)	32
1.1.	Objectives	32
1.2.	Equipment	32

1.3.	Procedure	33
1.4.	Practice questions	35
Practica	I 2. Dynamic Spirometry	37
2.1.	Objectives	37
2.2.	Equipment	37
2.3.	Procedure	38
2.4.	The Flow-volume loop (FVL)	39
2.5.	Volume-time graph (FEV₁ curve)	41
2.6.	Normal values	42
2.7.	Diagnostic Differences between Obstructive and Restrictive Airway	
Disea	ses	43
2.8.	Practice questions	44
2.9.	Further resources	46
2.10.	Summary	46
Reference	2	17

List of Abbreviations

CBC Complete blood count

ERV expiratory reserve volume

ESR Erythrocyte sedimentation

FEV₁ Forced expiratory volume at one second of forced vital

capacity

fl Femtoliter

FVC Forced vital capacity

FVL Flow-volume loop

Hb Hemoglobin

Hct Hematocrit

IC Inspiratory capacity

IRV Inspiratory reserve volume

MCH Mean corpuscular hemoglobin

MCHC Mean corpuscular hemoglobin concentration

MCV Mean corpuscular volume

MEF_{50%} Mean expiratory flow at 50% of forced vital capacity

PCV Packed cell volume

PEFR Peak expiratory flow rate

pg Picogram

PIFR Peak inspiratory flow rate

RBC Red blood cell

TLC Total lung capacity

TV Tidal volume

VC Vital capacity

WBC White blood cell

List of Figures

Figure 1. Procedure of drawing blood into capillary tubes. (A) The tip of a finger is
pricked using a lancet. (B) A heparinized capillary tube is gently brought in contact with
the blood drop forming on the tip of the finger14
Figure 2. Westergren's tubes and apparatus15
Figure 3. Drawing blood from a superficial vein in the antecubital fossa
Figure 4. Steps to doing a thin blood film. (A) a drop of anticoagulated blood is placed
near the edge of a microscope slide. The edge of a second slide is then placed at an
angle of 45° over the 1 st slide. (B) The edge of the second slide is slowly brought
towards the blood drop. Once it touched the blood drop, the blood drop is allowed to
spread along the edge of the 2 nd slide, then the 2 nd slide is pulled gently but swiftly over
the forst one to spread the blood over the 1 st slide creating a thin blood film (C)21
Figure 5. Schematic representation of the types of white blood cells that can be seen in
a peripheral blood film with a representative image of each WBC subset21
Figure 6. Blood group determination. (A) O ⁺ blood sample, (B) B ⁺ blood sample and (C)
A ⁺ blood sample24
Figure 7. Simple (volumetric) spirometer32
Figure 8. A spirogram recording. The deflection of the pen upwards or downwards with
each phase of respiration is dependent on machine mechanics and is subject to
variability
Figure 9. Simple lung volumes. The subject breaths through a mouthpiece while a nose
clip is placed on the nose to avoid air escaping through it. While breathing, air moves in
and out of the spirometer chamber causing displacement in the pen attached to it
surface. The moving pen draws the spirometry graph on the kymograph. The degree of
displacement is proportional to the volume of air moving in and out of the lungs. With
proper calibration, the volume of air moving in and out of the lungs can be calculated. 34
Figure 10. Automated spirometer38
Figure 11. Dynamic spirometry graphs. (A) Flow-volume loop, (B) FEV1 curve 39
Figure 12. A normal flow-volume loop40
Figure 13. Shows a normal compared to FVLs of patients with obstructive and restrictive
pulmonary disorders41
Figure 14. A normal volume-time graph (FEV ₁ curve)42
Figure 15. Typical spirograms and FVL in different ventilatory conditions (5),

List of Tables

Table 1. First year medical students' physiology practical sessions	9
Table 2. Faculty and staff members involved in teaching physiology practicals	10
Table 3. Normal FEV1 values (% predicted).	42
Table 4. Pulmonary function test interpretation.	43

Introduction

Welcome to your first year in medical school. This manual has been prepared as a

reference to help medical students navigate their way through 1st year physiology

practical sessions. The physiology practical sessions are part of the physiology

curriculum that is ingrained in the block system of the first two years of medical school.

Its aim is to provide a practical aspect to some of the physiological concepts learned

during mainstream lectures allowing students to have hands on experience that will

strengthen their understanding of the physiological concepts. Practical sessions will also

help them apply the knowledge learned in the classroom in a safe environment.

This manual is meant to be a guide for students providing the structure and topics

covered in physiology practical sessions. However, students are encouraged to look for

information and broaden their knowledge using other resources.

To make the best of the practical sessions, students are advised to attend the sessions

on time and prepare by reading related lecture material prior to the practical sessions.

During the sessions, students are encouraged to engage actively and take the

opportunity to get hands on experience whenever it is feasible.

Wish you all the best!

Physiology Practical Team

8

Overview of 1st Year Medical Student Physiology Practical Sessions

The table below provides an overview of the structure and contents of Physiology practical sessions during the first year in medical school by showing the number and title of the lab sessions provided in each block of the 1st year medical curriculum.

Table 1. First year medical students' physiology practical sessions.

Block name	Block	Number of sessions Session title	
	duration		
Foundation block	≈ 8 weeks	3 (2-hour-sessions)	CBC & ESR
			WBC & differential
			Bleeding and clotting times and
			blood groups
Musculoskeletal	≈ 5 weeks	None	None
block			
Respiratory block	≈ 6 weeks	2 (2-hour-sessions)	Lung volumes & capacities
			Dynamic spirometry
Cardiovascular	≈ 5 weeks	4 (2-hour-sessions)	The electrocardiogram (ECG)
block			Jugular venous and carotid
			arterial pressure recordings
			Heart sounds
			Measurement of arterial blood
			pressure
Renal block	≈ 5 weeks	4 (2-hour-sessions)	Glomerular filtration and renal
			clearance
			Diuresis-1
			Diuresis-2
			Acid-base balance

Faculty and Staff Members Involved in Physiology Practical Teaching

Table 2. Faculty and staff members involved in teaching physiology practicals.

Male members	Contact details	Female members	Contact details
Dr. Taj	athalepota@yahoo.com	Dr. Ola Hilmi	omawlana@ksu.edu.sa
Dr. Mustafa	mkmemon@gmail.com	Dr. Reem Al-Twairgi raltaweraqi@ksu.ed	
Dr. Yahya		Mrs. Sulafa Al-	Snaak-
Di. Tanya		Thubaiti	2009@hotmail.com
Mr. Jarouni		Mrs. Shurooq Al-Saidi	snaalsuod@gmail.com
Mr. Timhar		Mrs. Rahma	

Advice for Students

- Show up to physiology practical sessions on time.
- Read related lecture material prior to the laboratory session.
- Each laboratory session starts with a pre-lab lecture that serves to explain the objectives and procedure to be done in the lab. Listen attentively to these lectures.
- Engage actively in the laboratory activity and take every opportunity to get hands on experience whenever possible.
- Answer the sample questions provided at the end of each lesson.
- Do not depend solely of this guide for information.

Chapter 1: Foundation Block Physiology Practical

There are 3 practical physiology sessions during the foundation block which are all concerned with aspects related to blood physiology. The sessions are as follows;

- 1. Complete blood count (CBC) and red cell indices, erythrocyte sedimentation rate (ESR) and hematocrit (Hct).
- 2. White blood cell count (WBC) & differential.
- 3. Bleeding and clotting times and blood groups.

Practical 1. Complete Blood Count (CBC), Red blood cell indices, Erythrocyte Sedimentation Rate (ESR) and Hematocrit (Hct)

1.1. Objectives

At the end of session, the students should be familiar with:

- The procedures used for taking both capillary and venous blood.
- The methods used to measure the ESR and Hct.
- The normal values recorded when making these measurements.
- The method used to get CBC and assess red blood cell indices, including; mean corpuscular volume (MCV), mean corpuscular hemoglobin (MCH) and mean corpuscular hemoglobin concentration (MCHC).

1.2. Equipment

- Coulter analyzer.
- EDTA tubes.
- Lancets.
- Tourniquet.
- Alcohol swabs.
- Heparinized capillary tubes.
- Plasticine.
- A centrifuge.
- Micro-hematocrit reader.

1.3. Procedure

1.3.1. Measurement of Hct (or packed cell volume "PCV")

Blood is drawn from capillaries in order to measure Hct using micro-hematocrit reader.






Figure 1. Procedure of drawing blood into capillary tubes. (A) The tip of a finger is pricked using a lancet. (B) A heparinized capillary tube is gently brought in contact with the blood drop forming on the tip of the finger.

- 1. Clean the area of the skin of a finger-tip or an ear lobe with a sterilized alcohol swab.
- 2. Prick the skin using the pen lancet, Fig-1A.
- 3. Discard the first drop of blood, because it is mixed with tissue fluid.
- 4. Allow the second drop of blood to be formed and allow it to become large enough to fill 75% of the heparinized capillary tube by the capillary action when it is brought closer to the blood, Fig-1B. Apply only gentle pressure beneath the pricked skin to help the flow of blood, because if more pronounced pressure is exerted, blood is likely to be diluted with interstitial fluid.
- 5. Seal one end of the capillary tube with plasticine.
- 6. Repeat above steps 1 − 5 to collect several capillary blood samples.
- 7. Put all the capillary blood samples in a centrifuge machine for 5 minutes at the speed of 3000-4000 RPM to separate plasma from cells.
- Once centrifuged, take one of the capillary blood samples to see the cells have been packed at the bottom of the tube and the light-weight clear plasma visible above the cells.

9. Hct or PCV can then be determined as a percentage of the total volume using micro-hematocrit reader.

1.3.2. Measuring the erythrocyte sedimentation rate (ESR)

To measure ESR, the following equipment will be needed;

- Westergren's sedimentation apparatus, Fig-2.
- EDTA tubes.
- Disposable sterile syringes and needles.



Figure 2. Westergren's tubes and apparatus.

1.3.2.1. Procedure

- 1. Using a sterile syringe, draw 1.6 ml of blood from a suitable vein.
- 2. Transfer the blood to a test tube containing EDTA to prevent clotting.
- 3. Fill the Westergren's tube with blood up to the zero mark.
- 4. Place the tube upright in the stand and leave like this for one hour.

5. Note down the depth of the column of clear plasma at the top of red blood cells in the tube after one hour. This will be E.S.R. reading.

Normally the value of E.S.R. ranges from 0mm to 7 mm and it is slightly higher in females than males due to less number of red blood cells.

1.3.3. Counting peripheral blood cells using the Coulter analyzer

Blood is drawn from a superficial vein in the antecubital fossa using a needle attached to a syringe.




Figure 3. Drawing blood from a superficial vein in the antecubital fossa.

- 1. Clean the area of the skin to be pricked. Usually the blood is drawn from median cubital vein in front of the elbow joint to collect venous sample, Fig-3.
- 2. Apply the tourniquet above the elbow joint to impede the flow of venous blood towards the heart for a while.
- 3. Use a disposable syringe to draw the blood from the vein.
- 4. Immediately transfer the collected blood from the syringe to EDTA anticoagulated tube to prevent blood from clotting.

- 5. Activate the Coulter analyzer machine and a probe will move across and down into aspirate position. The aspiration syringe draws 12 μ I of whole blood into the probe.
- 6. The Coulter Analyzer makes the necessary dilutions with the reagents automatically and accurately counts and measures the sizes of cells by detecting and measuring changes in electrical resistance when a particle (such as cell) in the conductive liquid passes through a small aperture. As each cell goes through the aperture, it impedes the current and causes a measurable pulse. The number of pulses signals the number of particles. The height of each pulse is proportional to the volume of that cell.
- 7. Finally all the hematological values are reported and printed.

1.3.4. Calculation of red blood cell indices

1.3.4.1. Mean corpuscular volume (MCV)

This is the average volume of a red blood cell in an individual measured in *femtoliters* (fl). MCV can be calculated from the Hct and total red blood cell (RBC) count using the following formula:

$$MCV = Packed\ Cell\ Volume\ x \frac{10}{RBC\ Count}$$

Normally, MCV ranges between 78-98 fl. A low MCV denotes smaller than normal RBCs which are then called *microcytes*. Whereas, a high MCV denotes a larger than normal RBCs which are known as *macrocytes*.

1.3.4.2. Mean corpuscular hemoglobin (MCH)

This is the average weight of hemoglobin in a single red blood cell measured in picograms (pg).

$$MCH = Hemoglobin\ Concentration\ x \frac{10}{RBC\ Count}$$

Under normal conditions, the MCH ranges between 27–32 pg.

Low MCH denotes lower than normal hemoglobin weight in an RBC which is known as a hypochromic RBC. While, a high MCH denotes a higher than normal hemoglobin weight in an RBC which is called **hyperchromic** RBC.

1.3.4.3. Mean corpuscular hemoglobin concentration (MCHC)

This is the concentration of hemoglobin per 100 ml of red blood cell measured in grams/deciliters (g/dl).

$$\textit{MCHC} = \textit{Hemoglobin Concentration } x \frac{100}{\textit{Packed Cell Volume}}$$

Normally, the MCHC ranges between 32–36 g/dl. A MCHC value below normal suggests iron deficiency anemia.

1.4. Essential terminology

Find the meaning of the following medical terms;

- Polycythemia.
- Anemia.
- · Leucocytosis.
- Leucopenia.
- Thrombocytosis.
- Thrombocytopenia.

1.5	. P	rac	tice	au	esti	ons
		· uu		ЧЧ	00 11	00

1. What is the clinical importance of knowing the red blood cell indices?

2. Discuss briefly the etiological classification of anemia?

3. Peripheral blood parameters of two adult males (subject A and subject B) are shown in the table below. Using the information shown in the table answer the questions a & b.

Laboratory parameter	Subject A	Subject B
RBC count	3.6 X 10 ⁶ / mm ³	2.5 X 10 ⁶ / mm ³
Hb concentration	7.2 g/dl	8 g/dl
Hct	25%	25%

- a. Calculate MCV, MCH and MCHC for each of these subjects.
- b. What are the red blood cell abnormalities seen in these men. List possible causes for each of these abnormalities?
- 4. What is meant by rouleaux formation? And why does rapid rouleaux formation increase the E.S.R?
- 5. What is the clinical significance of E.S.R.?
- 6. What conditions are associated with an increased E.S.R.?

Practical 2. White Blood Cell Count (WBC) and Differential

2.1. Objectives

At the end of this session, students should be able to:

- Identify the different types of white blood cells under the microscope.
- Describe the normal values expected for each leucocyte subset.
- Understand the clinical relevance of the differential leucocyte count in disease diagnosis.

2.2. Equipment

- Light microscope with an oil immersion objective.
- Mineral or cedarwood oil.
- Wright's stain.
- Microscope slides.

2.3. Procedure

- 1. Venous blood is drawn into anticoagulant EDTA tubes.
- 2. After gentle mixing, a drop of blood is aspirated using a pipette and is placed at one end of a labelled slide. Make sure that the blood drop is as small as possible.
- 3. Then using another slide, place the second slide at a 45° angle over the first slide and slide it gently toward the blood drop.
- 4. Once the edge of the second slide touches the blood drop, allow the drop to spread along the edge of that slide then gently but swiftly pull the top slide over the bottom one to spread the blood drop over the bottom slide creating a thin film of blood, Fig-4.
- 5. Allow it to dry, and then stain the film using Wright's stain.
- 6. Set the stained blood film under the oil immersion objective in an electron microscope.

7. Identify various types of white blood cells according to their histological characteristics.




Figure 4. Steps to doing a thin blood film. (A) a drop of anticoagulated blood is placed near the edge of a microscope slide. The edge of a second slide is then placed at an angle of 45° over the 1st slide. (B) The edge of the second slide is slowly brought towards the blood drop. Once it touched the blood drop, the blood drop is allowed to spread along the edge of the 2nd slide, then the 2nd slide is pulled gently but swiftly over the forst one to spread the blood over the 1st slide creating a thin blood film (C).




Figure 5. Schematic representation of the types of white blood cells that can be seen in a peripheral blood film with a representative image of each WBC subset.

2.4. Practice questions

1. Fill the table below by writing the histologic features of each WBC subset, its normal value in the blood and the disease conditions in which they may be elevated.

Cell type	Neutrophil	Eosinophil	Basophil	Monocyte	Lymphocyte
Histologic					
features					
Normal					
value in					
blood					
Conditions					
causing an					
increase in					
its level					

2. What stains are used in the preparation of blood films?			

Practical 3. Blood groups, Bleeding & Clotting Times

3.1. Objectives

At the end of the session, the students should be able to:

- Understand and practice the method used in determining blood groups.
- Be familiar with the ABO and Rh systems of blood grouping and explain their importance in blood transfusion.
- Discuss the normal ranges of bleeding time and clotting time and determine their own values experimentally.
- Recognize the importance of bleeding time and clotting time in hemostasis.

3.2. Determination of blood groups

3.2.1. Equipment

- High titer anti-A, anti-B and anti-D sera.
- A microscope.
- Tooth picks.
- Microscope slides.
- Alcohol swabs.
- Lancet.

3.2.2. Procedure

- 1. Take 3 microscope slides and label them clearly as "A", "B" and "D".
- 2. Sterilize the fingertip with an alcohol swab.
- 3. Prick the finger using a lancet and place one drop of blood in each of the 3 microscope slides.

- 4. Quickly add a drop of anti-A, anti-B and anti-D sera to slides labeled as "A", "B" and "D" respectively.
- 5. Stir the mixture on each slide with the help of different pieces of tooth picks for a minute or two.
- 6. Examine the mixtures carefully for the signs of red blood cell agglutination. When red blood cells clump together (agglutination), they have a speckled or peppered appearance. If there is a doubt, examine the slides using the low power of a microscope. An example is shown in Fig-6.




Figure 6. Blood group determination. (A) O^+ blood sample, (B) B^+ blood sample and (C) A^+ blood sample.

3.2.3. Practice questions

1. What are the agglutinogens and agglutinins found in people with different blood groups in ABO system?

Blood group	Agglutinogen	Agglutinin
Α		
В		
AB		
0		

2. For each blood group in the table, write to which blood group/s can it donate blood and from which blood group/s can it receive blood in the setting of a blood transfusion.

Blood Group	Can give blood to	Can receive blood from
AB ⁺		
AB		
AB -		
A⁺		
Α -		
B⁺		
B ⁻		
O ⁺		
0.		

/hat is the distribution o	of the ABO and Rh blood groups in Saudi Arabi
Blood group	Percent in population
0+	
A+	
B+	
AB+	
0-	
A-	
B-	
AB-	
oes the distribution of b	olood groups in Saudi Arabia differ from that fo
rest of the world? /hat is hemolytic diseas	e of the newborn?

3.3. Determination of clotting time

3.3.1. Equipment

- Capillary tubes.
- Petri-dish.
- Alcohol swabs.
- Lancets.
- Plasticine.
- Water bath set at 37°C.
- A watch.

3.3.2. Procedure

- 1. Prick a finger of the subject observing the usual precautions and note the time at which the prick is made.
- 2. Wipe away the first drop of blood.
- 3. Then while the blood is still freely flowing, place one end of the capillary tube on it and let the tube fill with it by the capillary action.
- 4. Close both ends of this filled capillary tube with the plasticine.
- 5. Place this capillary tube in the water bath.
- 6. Repeat all the above steps with many capillary tubes.
- 7. Two minutes after making the prick, break a capillary tube and separate the two halves slowly and look for a thread like clot between the two broken halves of the tube.
- 8. Repeat step 7 at 30 seconds interval with the remaining tubes until you see a thread-like clot between the broken halves of one of the capillary tubes.
- 9. Note the time. The time from pricking the finger to the appearance of the clot is the clotting time.

3.3.	3. Practice questions
1.	What is the normal range for clotting time?
2.	What is/are the clinical condition/s in which the clotting time is greater than normal?
3.	Can you name some substances that are used as anti-coagulants.
4.	What is the clinical significance of clotting time?
5.	What is the source of heparin in the body?

3.4. Determination of bleeding time

3.4.1. Equipment

- Blotting paper.
- Stop watch.
- Alcohol swabs.
- Lancets.

3.4.2. Procedure

- 1. Prick a finger of the subject observing the usual precautions and note the time at which the prick is made. (The pricked skin should not be touched until the experiment is over.)
- 2. Apply a piece of filter paper (blotting paper) to the emerging drop of blood from the pricked skin every 30 seconds until the bleeding stops.
- 3. Note the time when the bleeding stops. The time from pricking the finger to the stop of bleeding is the bleeding time.

3.4.3. Practice questions

1.	What is the normal range of bleeding time?	
2.	An abnormality in which blood cell may prolong the bleeding time?	

3. Name one condition in which bleeding time is prolonged (increased)?

Chapter 2: Respiratory Block Physiology Practical

There are 2 practical physiology sessions during the respiratory block which include;

- 1. Static spirometry (lung volumes & capacities).
- 2. Dynamic spirometry.

Spirometry is one of the essential pulmonary function tests performed in clinical practice. It is concerned with the measurement of flow and volume of air entering and leaving the lungs. Two major types of spirometry measurements are usually performed: simple and dynamic. Simple spirometry is used for the determination of lung volumes and capacities, whereas dynamic spirometry measures the flow of air moving in and out of the lungs. In the following sections you will be introduced to the procedures used for measuring both simple and dynamic spirometry, their indications and result interpretation.

Practical 1. Simple Spirometry (Lung Volumes and Capacities)

1.1. Objectives

At the end of this session, students are expected to:

- 1. Describe how a bell-type spirometer is used to measure lung volumes and capacities.
- 2. List and define the different lung volumes and capacities.
- 3. State the normal values of each lung volume and capacity.
- 4. Discuss the physiological and pathological factors that may affect the different lung volumes and capacities.

1.2. Equipment

- 1. Simple spirometer (many types are available, Bell-type spirometer or water-gauge spirometer), Fig-7. It would be best if students acquaint themselves with the type used in the lab.
- 2. Nose clip.
- 3. Disposable mouth piece.



Figure 7. Simple (volumetric) spirometer.

1.3. Procedure

- 1. Insert the mouthpiece in the subject's mouth so that its edges lie between the subject's lips and gums.
- 2. Place the nose clip on the subject's nose to avoid air escaping through the nose.
- 3. Ask the subject to take normal breaths through the mouthpiece for a short while.
- 4. After recording few normal breaths, ask the subject to take a deep forceful inspiration filling their lungs to their maximum ability followed by gentle exhalation. After that, the subject can resume normal breathing.
- 5. After a few normal breaths, ask the subject to expire quickly, forcibly and as completely as possible. Once this forceful expiration is complete, the subject inhales and resumes normal breathing.
- 6. Finally, ask the subject to take a deep forceful inspiration followed immediately by a maximum, quick and forceful expiration. Once this is complete, ask the subject to breath normally for a short time.
- 7. The spirogram is recorded on a moving drum, Fig-8. An example of how the recording is done is shown in Fig-9.




Figure 8. A spirogram recording. The deflection of the pen upwards or downwards with each phase of respiration is dependent on machine mechanics and is subject to variability.




Figure 9. Simple lung volumes. The subject breaths through a mouthpiece while a nose clip is placed on the nose to avoid air escaping through it. While breathing, air moves in and out of the spirometer chamber causing displacement in the pen attached to it surface. The moving pen draws the spirometry graph on the kymograph. The degree of displacement is proportional to the volume of air moving in and out of the lungs. With proper calibration, the volume of air moving in and out of the lungs can be calculated.


Note to students

Depending on the mechanics of the machine used for simple spirometry measurements, the inspiratory/expiratory curves may be recorded upwards or downwards. The direction of inspiration and/or expiration will always be highlighted in any simple spirogram recording.

1.4.	Pra	ıctice	au	estio	ns
		.00.00	ЧЧ		

1.	Define the following terms and state/calculate their values from the data collected in the lab: a. Tidal volume (TV). b. Expiratory reserve volume (ERV). c. Inspiratory reserve volume (IRV). d. Vital capacity (VC). e. Inspiratory capacity (IC).
2.	Name a few physiological factors that may influence lung volumes and capacities and how do they exert their effects?
3.	Lung volumes and capacities are altered in a variety of pathological conditions. Name a few and explain how do these conditions bring about the changes are in lung volumes and capacities.
4.	What is the physiological significance of the residual volume and the functional residual capacity?
5.	Residual volume cannot be directly measured by spirometry. What is the technique that can be used to measure it? Explain how it works.

6. Using a simple ruler and the calibration provided in the graph, calculate the TV, IRV, ERV and VC from the graph below.



Parameter	Volume in liters
TV	
IRV	
ERV	
VC	

Practical 2. Dynamic Spirometry

2.1. Objectives

At the end of this session, students should be able to:

- 1. Perform a dynamic spirometry test on a fellow student.
- 2. Describe the two graphs recorded by dynamic spirometry, namely: flow-volume loop (FVL) and the volume-time curve (forced expiratory curve "FEV₁" curve).
- 3. Analyze the components of each graph; FVL and FEV₁ and describe the characteristics of a normal FVL and FEV₁ graphs.
- 4. Calculate the forced expiratory volume in the first second (FEV₁) and forced vital capacity (FVC) and the FEV₁/FVC ratio from the FEV₁ curve.
- 5. Calculate the FVC, peak expiratory flow rate (PEFR), peak inspiratory flow rate (PIFR) and maximal expiratory flow rate at 50% of the forced vital capacity (MEF₅₀).
- 6. Discuss the indications of dynamic spirometry in clinical practice.
- 7. State the normal values for FEV₁, FVC and the FEV₁/FVC ratio.
- 8. State the normal values of FVC, PEFR, PIFR and MEF₅₀ in FVL.
- 9. Describe the expected changes in FVL and FEV₁ curve in obstructive vs restrictive lung disease conditions.
- 10. Describe the expected changes in FEV₁, FVC and the FEV₁/FVC ratio in obstructive vs restrictive lung disease conditions.
- 11. Describe the expected changes in FVC, PEFR, PIFR and MEF₅₀ in obstructive vs restrictive lung disease conditions.

2.2. Equipment

- 1. Dynamic spirometer, Fig-10.
- 2. Nose clip.
- 3. Disposable mouth piece.




Figure 10. Automated spirometer.

2.3. Procedure

- 1. Insert a new disposable mouthpiece into the flow sensor (SP-250).
- 2. Hold the sensor in an upright position.
- 3. Insert the mouthpiece in the oral cavity (mouth) and seal the lips tightly around the mouthpiece.
- 4. Place the nose clip on the subject's nose to avoid air escaping through nostrils.
- 5. While subject is standing, allow him/her to breathe normally through mouthpiece, approximately 3 normal breaths to record TV.
- 6. Then ask the subject to inhale as deep as possible and then follow it with a fast and forceful expiration. The expiration should be as fast and forceful as possible and it should continue until the subject is unable to blow out anymore.
- 7. Two types of graphs may be recorded, Fig-11.




Figure 11. Dynamic spirometry graphs. (A) Flow-volume loop, (B) FEV1 curve.

2.4. The Flow-volume loop (FVL)

The FVL depicts the relationship between flow and volume under maximal effort of inspiration and expiration. The shape of the loop depends on the mechanical properties of the lung and may help in the diagnosis of ventilatory dysfunction. Fig-12 shows a normal FVL. The normal expiratory portion of a well-performed flow-volume loop is characterized by a rapid increase to the peak flow rate, followed by a nearly linear decrease in flow as the subject exhales toward residual volume. While normally the inspiratory portion shows a symmetric, saddle-shaped curve (1). The parameters that each student need to be familiar with and able to extrapolate from the FVL are: the peak expiratory flow rate (PEFR), peak inspiratory flow rate (PIFR), forced vital capacity (FVC) and maximum expiratory flow at the half-way point in the forced expiratory maneuver (MEF_{50%})-Fig-12.

Fig-13 shows FVL in normal compared to obstructive and restrictive pulmonary disorders.






Figure 12. A normal flow-volume loop. The expiratory portion of the loop lies above the x-axis while the inspiratory potion of the loop lies below the x-axis. PEFR represents the maximal flow achieved during forced expiration while the PIFR represents the maximal flow achieved during inspiration. The FVC is the total expiratory volume from a maximally forced expiration maneuver. $MEF_{50\%}$ is determined from the graph by first establishing the point at which 50% of the vital capacity has been expired (i.e. 2L in the graph above). A line perpendicular to the x-axis (volume axis) is drawn from this point towards the expiratory curve (dotted line no. 1). At the point of intersection between dotted line no. 1 and the expiratory curve, another line is drawn (dotted line no. 2) perpendicular to dotted line no. 1 towards the y-axis (flow axis). The point of intersection of the y-axis with dotted line 2 represents the $MEF_{50\%}$.








Figure 13. Shows a normal compared to FVLs of patients with obstructive and restrictive pulmonary disorders (reprinted from (2).

2.5. Volume-time graph (FEV₁ curve)

The volume-time graph or the FEV_1 curve depicts changes in volume (x-axis) against time (y-axis). Three main parameters are measured, namely, FVC, forced expiratory volume in the 1st second (FEV₁) and the ratio between these two numbers (FEV₁/FVC), Fig-14. When performing the test, one must ensure that the FEV₁ curve has reached a plateau and that expiration is maintained for at least 6 seconds (3, 4).




Figure 14. A normal volume-time graph (FEV₁ curve). The FVC represents the volume of air expired at the plateau. FEV_1 is the volume of air expired after 1 second of forced expiratory effort. Under normal conditions, more than 70% of the FVC is expired during the first second of expiration and this is what the ratio FEV_1/FVC reflects.

2.6. Normal values

Normal values are usually reported in 2 ways: as a volume measurement (ml or L of air), or as a percentage of the predicted normative or expected value for that patient's age, height, gender, and race from data obtained in the National Health and Nutrition Examination Survey III (NHANES III) (1).

Table 3. Normal FEV1 values (% predicted).

Parameter	Normal value (ATS/ERS)
FEV ₁	≥ 70% (% predicted FEV1)
FVC	≥ 70% (% predicted FVC)
FEV ₁ /FVC ratio	≥ 70% (0.7)

ATS=American Thoracic Society,

ERS=European Respiratory Society

2.7. Diagnostic Differences between Obstructive and Restrictive Airway Diseases

Using spirometry, pulmonary disorders may be categorized into:

- Obstructive.
- Restrictive.
- Mixed.

Obstructive pulmonary disorders are characterized by expiratory airflow limitation and can be seen as a disproportionate reduction in FEV_1 as compared to FVC. While restrictive pulmonary disorders are characterized by a reduction in FVC. Table-2 shows the characteristic findings in FEV_1 curve in the different ventilatory defects.

Table 4. Pulmonary function test interpretation.

	Obstructive pattern	Restrictive pattern	Mixed pattern		
FEV ₁	$\downarrow\downarrow\downarrow$	Normal or ↓	$\downarrow\downarrow$		
FVC	Normal or ↓	$\downarrow\downarrow\downarrow$	$\downarrow\downarrow$		
FEV ₁ /FVC (FEV ₁ %)	< 0.7 (70%)	Normal or > 0.7 (70%)	variable		

2.8. Practice questions

1. From the \mbox{FEV}_1 curve produced in the lab, what is the value of the following:

	Value					
Parameter	Litres	% predicted				
FEV ₁						
FVC						
FEV ₁ /FVC ratio						

2	. What is the expected normal value for FEV₁ in a normal person?
3	. How long does it take for healthy subjects to expire approximately 70% of their vital capacity?
4	Briefly explain what happens to FVC, FEV_1 and FEV_1 % measurements in patients with obstructive and restrictive lung diseases.

5.	From	the	flow	volume	loop	recorded,	what	is th	e value	of	the	follow	ing
	paran	neter	s:										

Parameter	Value	
	Litres	% predicted
PEFR		
PIFR		
FVC		
MEF _{50%}		

6	Briefly describe the important characteristics of the flow-volume curve recorded in a normal healthy person.
7.	Why is the force-independent part of the expiratory loop curvilinear in obstructive lung disease?
8.	What is the clinical significance of MEF50 measurements?

2.9. Further resources

- Paraskeva et al. Spirometry. 2011. Australian Family Physician. 40 (4): 216-219.
- Johnson et al. A stepwise approach to the interpretation of pulmonary function tests. 2014. American Family Physician. 89 (5): 359-366.

2.10. Summary

Fig-15 summarizes the findings seen in FVL and FEV_1 curve in abnormal ventilatory conditions compared to normal.




Figure 15. Typical spirograms and FVL in different ventilatory conditions (5).

References

- 1. Parker MJ. Interpreting spirometry: the basics. Otolaryngol Clin North Am. 2014;47(1):39-53.
- 2. Paraskeva MA, Borg BM, Naughton MT. Spirometry. Aust Fam Physician. 2011;40(4):216-219.
- 3. Barriero TJ, Rerillo I. An approach to interpreting spirometry. Am Fam Physician. 2004;69(5):1107-1114.
- 4. Johnson JD, Theurer WM. A stepwise approach to the interpretation of pulmonary function tests. Am Fam Physician. 2014;89(5):359-366.
- 5. Pierce R. Spirometry: an essential clinical measurement. Aust Fam Physician. 2005;34(7):535-539.