



# **DNA Extraction and Purification**

# Lab Equipment

Automatic pipettes



Vortex




Microcentrifuge



Water bath



UV-spectrophotometer




# Continued ...



Eppendorf tube




Cuvettes



Tips



Rack- test tube



Rack- eppendorf tube



# DNA Extraction

## Principle:

1. **Lysis of nucleated cells.**
2. **Removal of contaminants:** Any substance other than DNA, e.g., proteins.
3. **Measurements:** UV absorbance at 260nm and 280nm

**Purity of DNA solution:** 260/280 ratio


**DNA concentration:** Absorbance at 260nm



## *Continued ...*

### Steps:

1. **Lysis** of nucleated cells using lysis buffer.
2. **Binding** of DNA to the membrane of spin column.
3. **Wash:** using wash buffer.
4. **Elution** of pure DNA.







# **Spin Protocol of DNA Purification from Blood**

Video



<https://www.youtube.com/watch?v=gmNw6CWtN5k>



1. Pipette **20 $\mu$ l** protease.
2. Add **200 $\mu$ l** sample.
3. Add **200 $\mu$ l** Buffer **AL**.




4. Mix by pulse-vortex for 15s and incubate at **56°C** for **10 minutes**.
5. Briefly centrifuge.
6. Add **200 $\mu$ l** 96-100% ethanol and mix by pulse-vortex for 15s.
7. Briefly centrifuge.



8. Apply the mixture to the Mini spin column.
9. Centrifuge at **8000 rpm** for **1 minute**.



10. Discard the collection tube.
11. Place the column in a new 2ml collection tube.
12. Add **500 $\mu$ l** buffer **AW1**.
13. Centrifuge at **8000 rpm** for **1 minute**.





14. Place the column in a new 2ml collection tube.
15. Add **500 $\mu$ l** buffer **AW2**.
16. Centrifuge at **14000 rpm** for **3 minutes**.
17. Place the column in a new 2ml collection tube.
18. Centrifuge at **14000 rpm** for **1 minute**.




19. Place the column in a new 1.5ml tube.
20. Add **200 $\mu$ l** buffer **AE**.
21. Incubate at room temperature for **1 minute**.
22. Centrifuge at **8000 rpm** for **1 minute**.



**Add 400 $\mu$ l buffer AE to dilute the eluted DNA (3X dilution)**



23. Quantify the DNA concentration.







# **Quantification of the purified DNA**

# Measurements

Done by  
the spectrophotometer

- Measure the Absorbance at **260nm** and at **280nm**.
- **Assess the DNA purity:** 260/280 ratio  
(Accepted ratio: **1.7 - 1.9**)
- **Calculate DNA Conc.:** Provided  $A_{260} = 1.0$ , DNA is  $50 \mu\text{g/ml}$ , unknown DNA Conc. can be calculated by cross multiplication

$$A_{260} = 1.0$$

$$\text{DNA conc.} = 50 \mu\text{g/ml}$$

$$A_{260} = 0.5$$

$$\text{DNA conc. ?}$$

**Note:** In case of diluting the eluted sample, multiplies the final concentration by the dilution factor. This can be adjusted by the spectrophotometer.

# DNA Yield

**DNA Yield = DNA Volume x final DNA Conc.**


**Example:**

If you have

**Volume of DNA solution:** 200 $\mu$ l (0.2 ml)

**Final DNA Conc.:** 30  $\mu$ g/ml

**Then, the yield ( $\mu$ g) = 0.2 ml x 30  $\mu$ g/ml  
= 6.0  $\mu$ g**






# **Molecular Techniques and Applications**

***Note:***


*Almost all molecular biology techniques can be utilized for diagnosis and research*





# a. Amplification techniques:

## e.g. Polymerase Chain Reaction (PCR)




- 1 **Denaturation** at 94-96°C
- 2 **Annealing** at ~68°C
- 3 **Elongation** at ca. 72 °C



# ***Continued ...***

## ***Applications of PCR:***

- *Comparison of a normal gene with a mutant form of the gene.*
- *Detection of low-abundance nucleic acid sequences.*
- *Forensic analysis of DNA samples.*
- *Prenatal diagnosis.*



# Other examples of molecular techniques:

1. Restriction Fragment length polymorphism (RFLP).



<https://www.ncbi.nlm.nih.gov/probe/docs/techrflp/>

2. Southern blotting.



<http://www.onlinebiologynotes.com/southern-blotting-principle-procedure-application/>

