

King Saud University
College of Medicine
Clinical Biochemistry unit, Department of Pathology

Biochemistry Practical Class

Glucose Estimation In Blood & Urine

Endocrine System Block (2nd Year)

Multipurpose Lab | College of Medicine | Level 2



Endocrine Block

Biochemistry Practical

Practical biochemistry sessions are designed to teach the fundamental techniques used in biochemistry to medical students. In this session the students will have hands-on training in estimating glucose in blood and urine. This will enable them to understand how biochemical data is collected, analyzed and interpreted.

Plan of work:

- 1. Handouts will be distributed and a presentation on glucose estimation in blood and urine will be given for 15 minutes.
- 2. Hands-on practical session will be conducted following the presentation.

Estimation of Glucose in Blood

Principle of the method

Blood glucose is detected by a series of enzymatic reactions (oxidase and peroxidase) that ultimately form a colored product. The intensity of the color is directly proportional to the amount of glucose in blood. The color intensity is determined by measuring the absorbance at a wavelength of 546 nm using a spectrophotometer.

Procedure:

	Test	Standard	Blank
Reagent	2.5 ml	2.5 ml	2.5 ml
Blood sample	25 μl	-	
Standard	-	25 μl	-
H ₂ O	-	-	25 µl

- Mix and incubate for 10 min. at room temperature.
- Measure the absorbance (A) of the test and standard against blank at 546 nm.

Calculation:

Glucose conc. (mmol/l) =
$$\frac{A \text{ Test}}{A \text{ Standard}}$$
 X conc. of standard. (5.6 mmol/l)

Normal reference values for blood glucose:

3.9-6.1 mmol/L (70-110 mg/dL)

Note: To convert mmol/l to mg/dL multiply the results by 18

Urine analysis (using dipstick):

Principle:

- Dipsticks are plastic strips impregnated with chemical reagents that react with specific substances in the urine to produce color-coded visual results.
- They provide quick determination of pH, protein, glucose and ketones.
- The depth of color produced is directly proportional to the concentration of the substance in urine.
- Color controls are provided against which the actual color produced by the urine sample can be compared. The reaction times of the impregnated chemicals are standardized.

Procedure:

- Dip the dipstick in the urine sample provided and remove it immediately.
- Wipe off the excess urine.
- Read the color produced within 60 sec. by comparing with control chart.

Item	Observation	Comment
Protein	Pasitive 566	Proteinuvea
Glucose	Positive	hyrerglycemia
Ketones	Positive +++	Keto Acidosis
pН	Positive 6.0	6.0

Criteria for the diagnosis of diabetes

 A1C ≥6.5 percent. The test should be performed in a laboratory using a method that is NGSP certified and standardized to the DCCT assay.*

OR

2. FPG ≥126 mg/dL (7.0 mmol/L). Fasting is defined as no caloric intake for at least 8 h.*

OR

3. Two-hour plasma glucose ≥200 mg/dL (11.1 mmol/L) during an OGTT. The test should be performed as described by the World Health Organization, using a glucose load containing the equivalent of 75 g anhydrous glucose dissolved in water.*

OR

 In a patient with classic symptoms of hyperglycemia or hyperglycemic crisis, a random plasma glucose ≥200 mg/dL (11.1 mmol/L).

A1C: glycated hemoglobin; NGSP: National glycohemoglobin standardization program; DCCT: Diabetes control and complications trial; FPG: fasting plasma glucose; OGTT: oral glucose tolerance test.

* In the absence of unequivocal hyperglycemia, criteria 1-3 should be confirmed by repeat testing.

Reprinted with permission from: American Diabetes Association.

Standards of Medical Care in Diabetes Care 2010; 33:S11.

Copyright © 2010 American Diabetes Association.

Categories of increased risk for diabetes*

FPG 100-125 mg/dL (5.6-6.9 mmol/L) [IFG]

2-h PG on the 75-g OGTT 140-199 mg/dL (7.8-11.0 mmol/L) [IGT]

A1C 5.7-6.4 percent

FPG: fasting plasma glucose; IFG: impaired fasting glucose; PG: post glucose; OGTT: oral glucose tolerance test; IGT: impaired glucose tolerance; A1C: glycated hemoglobin.